

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Det Johnson #531 AGENCY LPESO
DATE 3-13-16 TIME 1431
INTOXILYZER SERIAL # 80-001031 LOCATION SALOME

- 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- 4. Air blank completed.
- 5. Calibration check completed. Test results 0. 100 AC.
- 6. Air blank completed.
- 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE _____
DPS Form Exh G-2 (Rev 05-01)

INTOXILYZER 8000
Location: SALOME
Serial Number: 80-001031
Core Version: 8105.48
03/13/2016 14:31:36

Standard Lot#: 08915100A2
Last Changed By: J. GOMEZ #220

QAS: M. JOHNSON #531
LA PAZ CSO

Operator: M. JOHNSON #531
LA PAZ CSO

Test	g/210L	Time
Air Blank	0.000	14:32:06
0.100 Cal Check	0.100	14:32:27
Air Blank	0.000	14:32:57

Preventative Maintenance
Performed -- Timer Reset