

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Dep Johnson #531 AGENCY LPCSO
DATE 2-13-16 TIME 1321
INTOXILYZER SERIAL # 80-001196 LOCATION EHRENBURG

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- ~~2.~~ 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
~~3.~~ 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
~~4.~~ 4. Air blank completed.
~~5.~~ 5. Calibration check completed. Test results 0. 096 AC.
~~6.~~ 6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE _____

INTOXILYZER 8000

Location: EHRENBURG SUBSTATION

Serial Number: 80-001196

Core Version: 8105.48

03/13/2016

13:23:50

Standard Lot#: 08915100A2

Last Changed By: E. SCHNEIDER #431

QAS: M. JOHNSON #531

LA PAZ CSO

Operator: M. JOHNSON #531

LA PAZ CSO

Test	g/210L	Time
Air Blank	0.000	13:24:25
0.100 Cal Check	0.096	13:24:46
Air Blank	0.000	13:25:16

Preventative Maintenance
Performed -- Timer Reset