

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST E. Schneider AGENCY LPCS0
DATE 10/15/15 TIME 2358
INTOXILYZER SERIAL # 80-001196 LOCATION Ehrenberg Substation

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .10d AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 096 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

INTOXILYZER 8000

Location: EHRENBURG SUBSTATION

Serial Number: 80-001196

Core Version: 8105.48

10/15/2015

23:58:49

Standard Lot#: 08915100A2

Last Changed By: E. SCHNEIDER #431

QAS: E. SCHNEIDER #431

LA PAZ CSO

Operator: E. SCHNEIDER #431

LA PAZ CSO

Test	g/210L	Time
Air Blank	0.000	23:59:34
0.100 Cal Check	0.096	23:59:54
Air Blank	0.000	00:00:24

Preventative Maintenance
Performed -- Timer Reset