

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST SGT. R. EPPS AGENCY LAPAZ SO

DATE 10/3/12 TIME 1047

INTOXILYZER SERIAL # 80-001198 LOCATION EMLENBURG

- (initials)* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (initials)* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (initials)* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (initials)* 4. Air blank completed.
- (initials)* 5. Calibration check completed. Test results 0. 099 AC.
- (initials)* 6. Air blank completed.
- (initials)* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE A. Lynn

INTOXILYZER 8000

Location: EHRENBERG STATION

Serial Number: 80-001198

Core Version: 8105.48

11/03/2012

10:47:36

Standard Lot#: 15712100A1

Last Changed By: R. EPPS #90

QAS: R. EPPS #90

LA PAZ CSO

Operator: R. EPPS #90

LA PAZ CSO

Test g/210L Time

Air Blank 0.000 10:48:17
0.100 Cal Check 0.099 10:48:37
Air Blank 0.000 10:49:07

Preventative Maintenance
Performed -- Timer Reset