

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)
ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST DEPUTY JOHN M. GOMEZ #220 AGENCY LA PAZ COUNTY SHERIFF'S OFFICE

DATE 8/22/2015 TIME 17:27 Hrs

INTOXILYZER SERIAL # 80 - 001200 LOCATION QPD 103/LPCSO BACK UP

(✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.

OR

Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34 C 0.2 C

(✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

(✓) 4. Air blank completed.

(✓) 5. Calibration check completed. Test results 0. 0.91 AC.

(✓) 6. Air blank completed.

(✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE _____

INTOXILYZER 8000.

Location: QTZ 103/LPCS BACK UP

Serial Number: 80-001200

Core Version: 8105.48

08/22/2015

17:27:35

Standard Lot#: 08915100A2

Last Changed By: R. ROSADO #557

QAS: J. GOMEZ #220

LA PAZ CSO

Operator: J. GOMEZ #220

LA PAZ CSO

Test	g/210L	Time
Air Blank	0.000	17:28:11
0.100 Cal Check	0.091	17:28:32
Air Blank	0.000	17:29:01

Preventative Maintenance
Performed -- Timer Reset

Handwritten signature and initials, possibly "R. Rosado", in black ink.