

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST DEPUTY JOHN M. GOMEZ #220 AGENCY LA PAZ COUNTY SHERIFF'S OFFICE

DATE 12 / 07 / 2011 TIME 22:12 Hrs

INTOXILYZER SERIAL # 80 - 003861 LOCATION LPCSO BOATHOUSE

(✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.

OR

Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34 C 0.2 C

(✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.


(✓) 4. Air blank completed.

(✓) 5. Calibration check completed. Test results 0.102 AC.

(✓) 6. Air blank completed.

(✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE \_\_\_\_\_



INTOXILYZER 8000

Location: LPCSO BOATHOUSE

Serial Number: 80-003861

Core Version: 8105.48

12/07/2011

22:12:12

-----  
Standard Lot#: 769860

Last Changed By: R. EPPS #90  
-----

QAS: J. GOMEZ #220

LA PAZ CSO  
-----

Operator: J. GOMEZ #220

LA PAZ CSO  
-----

Test	g/210L	Time
Air Blank	0.000	22:12:44
0.100 Cal Check	0.102	22:13:04
Air Blank	0.000	22:13:34

Preventative Maintenance  
Performed -- Timer Reset



**EXHIBIT G-4**  
**THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)**  
**ARIZONA DEPARTMENT OF PUBLIC SAFETY**  
**STANDARD QUALITY ASSURANCE PROCEDURES**  
**INTOXILYZER MODEL 8000**  
**STANDARD CALIBRATION CHECK PROCEDURE**

**QA SPECIALIST** DEPUTY JOHN M. GOMEZ #220      **AGENCY** LA PAZ COUNTY SHERIFF'S OFFICE

**DATE** 12 / 07 / 2011      **TIME** 22:15 Hrs

**INTOXILYZER SERIAL #** 80 - 003861      **LOCATION** LPCSO BOATHOUSE

(✓) 1. Display reads "PUSH BUTTON TO START"

**DIAGNOSTIC TESTS**

- (✓) 1. Clock time check.
- (✓) 2. Date check.

**OPERATIONAL TESTS**

- (✓) 1. Alcohol-free subject test result 0.000 AC.
- (✓) 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed
- (✓) 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.
- (✓) 4. Standard Calibration Check standard 0.100 AC Result 0.103 AC.

Instrument is operating properly and accurately. YES  NO

**COMMENTS**


---

---

---

---

**SIGNATURE** \_\_\_\_\_



INTOXILYZER 8000  
Location: LPCSO BOATHOUSE  
Serial Number: 80-003861  
Core Version: 8105.48  
12/07/2011

22:15:06

---

\* Standard Lot#: 769860  
Last Changed By: R. EPPS #90

---

QAS: J. GOMEZ #220  
LA PAZ CSO

---

Operator: J. GOMEZ #220  
LA PAZ CSO

---

Subject: SQAP, TEST  
DOB: 01/02/1975  
Sex: M Weight: 175

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	22:15:57
Diagnostic Test	Pass	22:16:32
Air Blank	0.000	22:16:59
0.100 Cal Check	0.103	22:17:20
Air Blank	0.000	22:17:50
Subject Test	0.000	22:18:10
Air Blank	ABT*	22:18:17
Air Blank	IPS**	22:18:24

\*Sequence Aborted

\*\*Improper Sample *July 20*

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: LPCSO BOATHOUSE  
Serial Number: 80-003861  
Core Version: 8105.48  
12/07/2011

22:18:44

Standard Lot#: 769860  
Last Changed By: R. EPPS #90

QAS: J. GOMEZ #220  
LA PAZ CSO

Operator: J. GOMEZ #220  
LA PAZ CSO

Subject: SQAP,TEST  
DOB: 01/02/1975  
Sex: M Weight: 175

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	22:19:37
Diagnostic Test	Pass	22:20:12
Air Blank	0.000	22:20:39
0.100 Cal Check	0.103	22:21:00
Air Blank	0.000	22:21:29
Subject Test	0.077*	22:24:37
Air Blank	0.000	22:25:08
Air Blank	0.000	22:25:36
Subject Test	XXX**	22:26:04
Air Blank	RFI***	22:26:12

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

\*\*\*RFI Detect

Not a Successfully Completed Test Sequence

**EXHIBIT G-2**  
**THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R-13-10-104 (A)**  
**ARIZONA DEPARTMENT OF PUBLIC SAFETY**  
**STANDARD QUALITY ASSURANCE PROCEDURES**  
**INTOXILYZER MODEL 8000**  
**STANDARD CALIBRATION CHECK PROCEDURE**

**QA SPECIALIST** DEPUTY JOHN M. GOMEZ #220      **AGENCY** LA PAZ COUNTY SHERIFF'S OFFICE

**DATE** 12 / 07 / 2011      **TIME** 22:26 Hrs

**INTOXILYZER SERIAL #** 80 - 003861      **LOCATION** LPCSO BOATHOUSE

(✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.

**OR**

Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34 C 0.2 C

(✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"

(✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.

(✓) 4. Air blank completed.

(✓) 5. Calibration check completed. Test results 0.102 AC.

(✓) 6. Air blank completed.

(✓) 7. Remove printed record. Attach the record to the completed checklist.

**SIGNATURE** \_\_\_\_\_



INTOXILYZER 8000  
Location: LPCSO BOATHOUSE  
Serial Number: 80-003861  
Core Version: 8105.48  
12/07/2011

22:26:59

-----  
Standard Lot#: 769860  
Last Changed By: R. EPPS #90  
-----

QAS: J. GOMEZ #220  
LA PAZ CSO  
-----

Operator: J. GOMEZ #220  
LA PAZ CSO  
-----

Test	g/210L	Time
Air Blank	0.000	22:27:29
0.100 Cal Check	0.102	22:27:50
Air Blank	0.000	22:28:20

Preventative Maintenance  
Performed -- Timer Reset

*J. Gomez*  
#220